# Blom-Singer®

voice restoration systems

## Blom-Singer® Laryngectomy Tube

USER Instructions For Use

**INHEALTH®**  
TECHNOLOGIES

**R1**

37453-01E

37453-01E Effective March 2010 / Gültig ab März 2010 / Efectivo a partir de marzo de 2010 / Date d'entrée en vigueur : mars 2010 / Valido da marzo 2010 / Met ingang van maart 2010 / Vigente em de março de 2010

Blom-Singer is a registered trademark of Hansa Medical Products, Inc.  
InHealth Technologies is a registered trademark of Helix Medical, LLC  
© 2010 InHealth Technologies, A Division of Helix Medical, LLC

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## **BLOM-SINGER® LARYNGECTOMY TUBE**

The Blom-Singer Laryngectomy Tube (diagram 1) is a durable, lightweight, silicone tube that is used to keep the tracheostoma open after laryngectomy. The Laryngectomy Tube is available in a variety of sizes and flexibly contours to the postoperative tracheal anatomy. Along the posterior (backside) it features a recessed area with smooth, rounded edges that can be modified (fenestrated) for use with a Blom-Singer Voice Prosthesis. A pre-fenestrated series (diagram 2) is also available.

**This product is intended for single patient use only. Not for use between patients.**

**Do not use the device if it becomes contaminated.**

### **INDICATIONS**

The Blom-Singer Laryngectomy Tube is designed to maintain tracheostoma patency after laryngectomy to avoid tracheostomal closure.

### **CONTRAINDICATIONS**

The Blom-Singer Laryngectomy Tube is a medical product that should be used only by professionals with experience and training in its use and care. Prolonged infection, especially when resistant to therapy, may contraindicate use of this tube. Radiation therapy and resultant inflammation, bleeding, discharge, or coughing is an additional contraindication.

### **HOW SUPPLIED**

The Blom-Singer Laryngectomy Tube is available in three product offerings: **Sterile** and **Nonsterile** Laryngectomy Tube (non-fenestrated), and a **Nonsterile Fenestrated** Laryngectomy Tube. A laryngectomy tube holder neck strap is available separately.

### **INSTRUCTIONS FOR USE**

The user's hands and the laryngectomy tube should always be thoroughly cleaned to help avoid introducing contaminants into the tracheostoma and respiratory tract. Removal and reinsertion

of the device should only be done with the user positioned in front of a mirror with a bright light focused directly on the stoma.

*Please refer to the diagrams at the back of this instruction manual.*

*The following procedure is recommended by Eric D. Blom, Ph.D.*

**1.** Select the largest diameter and length laryngectomy tube that is comfortably tolerated by the patient and permits easy breathing. Insert the smaller tip end of the tube completely into the tracheostoma until its larger funnel shaped end rests against the entrance to the tracheostoma.

**2.** A very small amount of medical grade, water-soluble lubricant applied to the rim of the tracheostoma may ease insertion of the laryngectomy tube. In some instances it may be necessary to slightly fold or compress the tip portion of the laryngectomy tube, forming a taper, to facilitate introduction into the tracheostoma. The tube is curved and should conform comfortably to the post-laryngectomy contour of the trachea.

**3.** The laryngectomy tube may be secured around the neck using a laryngectomy tube holder (sold separately). The straps are inserted into the openings in the side of the laryngectomy tube so that the tabs on the straps point away from the tube. The straps are passed around the back of the neck and adjusted so that the straps are secure. The straps should be loose enough to permit one finger breadth between the backside of the straps and the surface of the neck. After adjusting the straps to a comfortable tension, the user may choose to cut off some of the excess strap length.

## FENESTRATING INSTRUCTIONS

A Blom-Singer Laryngectomy Tube used simultaneously with a voice prosthesis usually will require fenestration of the posterior wall to facilitate unrestricted airflow for tracheoesophageal voice. A pre-fenestrated series is available for this purpose. On non-pre-fenestrated laryngectomy tubes, a recessed area has been molded into the laryngectomy tube in this area to facilitate secondary modification by a physician or other experienced and trained clinician. The tube should not be in place during the fenestration procedure. Do not make the fenestration with razors or scissors. Tearing or other damage to the device can result if

cutting is done improperly. The laryngectomy tube should **never** be structurally modified by the user (patient).

*The following steps are recommended by Eric D. Blom, Ph.D., for fenestrating a Blom-Singer Laryngectomy Tube:*

**1.** Standard aseptic techniques should be observed when handling the laryngectomy tube. With the voice prosthesis positioned correctly in the tracheoesophageal puncture, carefully insert the appropriate diameter and length laryngectomy tube following the procedures recommended in the Instructions for Use in this data sheet.

**2.** With the laryngectomy tube fully inserted, visualize through the transparent posterior wall of the laryngectomy tube the point of airflow entry to the voice prosthesis. In most patients this will be within the recessed area on the tube. Using a pen, carefully mark that point against the interior surface of the laryngectomy tube.

**3.** Carefully remove the laryngectomy tube from the tracheostoma. Using a Fenestration Punch that has been autoclave sterilized or thoroughly disinfected with an antiseptic prior to use, make a hole corresponding to the ink mark. Only a 1/4" punch should be used. The punch is positioned by sliding it down over the funnel-shaped end of the laryngectomy tube. After the fenestration is placed, remove the punched disc completely. Be certain that all particles of silicone cut by the punch are removed from the tube.

**4.** Thoroughly rinse the laryngectomy tube with water and wash with an antiseptic. It is recommended that the laryngectomy tube be reinserted wet. **Never** use facial or toilet tissue to dry a silicone laryngectomy tube as particles of lint or fabric may be deposited onto the laryngectomy tube that can be inhaled into the airway upon reinsertion. If a laryngectomy tube must be dried, use only a lint-free material and rinse the laryngectomy tube again in water before reinsertion. Carefully inspect for loose particles and structural integrity. **Never** use a laryngectomy tube that has tears or other potentially hazardous structural weaknesses.

**5.** Reinsert the clean, disinfected laryngectomy tube into the tracheostoma according to the Instructions for Use. Observe for

---

adequate fenestration-to-voice prosthesis alignment, and assess the patient's capacity to phonate. If alignment is inadequate, follow the fenestration procedure previously described and enlarge the size of the punch hole as required. Clean and inspect again as described above.

## PATIENT CLEANING INSTRUCTIONS

The user's hands should be thoroughly washed with an antiseptic before handling the laryngectomy tube. Remove the laryngectomy tube from neck and rinse it clean with tap water. Wash the laryngectomy tube with an antiseptic and rinse again.

Prior to reinsertion, inspect the laryngectomy tube for structural damage, such as cracks or tears, resulting from prolonged use. Do not use a tube that is damaged.

It is recommended that the laryngectomy tube be reinserted wet. **Never** use facial or toilet tissue to dry a silicone laryngectomy tube as particles of lint or fabric may be deposited onto the laryngectomy tube that can be inhaled into the airway upon reinsertion. If a laryngectomy tube must be dried, use only a lint-free material and rinse the laryngectomy tube again in water before reinsertion.

Some users experience deposits of bacteria or yeast (*Candida albicans*) on the surface of the laryngectomy tube indicating the need for more effective device disinfection. These deposits appear as patches. In order to prevent these yeast deposits, and to routinely disinfect the laryngectomy tube, alternately wearing and soaking two devices is recommended.

After removing and cleaning the tube as described above, place it in a small container of fresh 3% hydrogen peroxide overnight. After soaking, the tube should be thoroughly rinsed again with tap water, air dried completely, and then placed in a clean resealable plastic bag until its next use. **Note:** Soaking in hydrogen peroxide will not remove yeast that is already present on the laryngectomy tube. It is a preventive measure only. If yeast is present, the laryngectomy tube should be replaced.

**Never insert any cleaning device or punch into the laryngectomy tube while the tube is positioned in the neck.** Always clean the device *after* it has been removed.

Do not use solvents or petroleum-based products for cleaning or lubricating the device. These materials may damage silicone.

## **WARNINGS AND PRECAUTIONS**

Never use facial or toilet tissue to dry a silicone laryngectomy tube as particles of lint or fabric may be deposited onto the laryngectomy tube that can be inhaled into the airway upon reinsertion.

Never insert any cleaning device or punch into the laryngectomy tube while the tube is positioned in the neck. Always clean the device after it has been removed.

Never use a laryngectomy tube that has tears or other potentially hazardous structural weaknesses.

The tube should not be in place during the fenestration procedure. Do not make the fenestration with razors or scissors. Tearing or other damage to the device can result if cutting is done improperly.

Tube modifications should not be attempted by the patient.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician or licensed speech pathologist.

## **COMPLICATIONS**

The following complications have been identified to occur with the use of this device, but are infrequent. They include: Stoma contamination or infection, which may require removal of the laryngectomy tube and/or appropriate antibiotics; Tissue granulation around the stoma; Accidental inhalation of an inappropriately sized laryngectomy tube into the airway, which may require removal by a physician; Obstruction of the tube by mucus necessitating removal and cleaning; Structural failure from excessive usage and/or cleaning.

---

## **ORDERING INFORMATION**

### **USA**

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: [order@inhealth.com](mailto:order@inhealth.com) ORDER ON-LINE: [www.inhealth.com](http://www.inhealth.com) POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

### **International**

Please contact our customer service department for a distributor referral.

### **Consumer Affairs**

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: [productcomplaints@inhealth.com](mailto:productcomplaints@inhealth.com)

### **RETURNED GOODS POLICY**

All returned merchandise must have a Return Merchandise Authorization number (RMA), and must be unopened and undamaged. RMA numbers are obtained from InHealth Customer Service. Returns without an RMA will not be accepted. Products must be returned in unopened packages, with manufacturer's tamper evident seals intact to be accepted for replacement or credit. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 14 days. Special order products are nonreturnable. There is a 20% restocking fee on all returned merchandise.

### **LIMITED WARRANTY, LIMITED REMEDY AND LIMITATION OF LIABILITY**

InHealth Technologies, a Division of Helix Medical, LLC, warrants that the Product will be free from defects in workmanship and materials on the date of purchase. THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES OR REPRESENTATIONS, EXPRESS OR IMPLIED, ARISING BY LAW OR CUSTOM, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR

PURPOSE, SAID IMPLIED WARRANTIES BEING HEREBY EXPRESSLY DISCLAIMED. Purchaser's sole and exclusive remedy for a valid warranty claim is replacement of the Product. This remedy does not include the cost of installation, removal, dismantling, delivery or reinstallation. InHealth shall have no liability to the extent Products are or have been modified or subject to misuse, abuse or improper storage, installation or maintenance or damaged in transit. No returns will be accepted without prior approval and authorization. INHEALTH SHALL NOT BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT, CONSEQUENTIAL, PUNITIVE, OR EXEMPLARY DAMAGES OF ANY KIND.

## **GLOSSARY**

### **Airway**

The air passage by which air reaches the lungs.

### **Candida albicans**

A yeast fungus that is a commonly occurring condition that can infect the mouth or the skin and may multiply and overgrow on to the voice prosthesis.

### **Clinician**

A health professional, such as a physician, licensed speech-language pathologist, or nurse who is involved in clinical practice.

### **Contraindications**

A reason that makes it inadvisable to prescribe a particular procedure or treatment.

### **Esophagus**

The muscular membranous tube for the passage of food from the pharynx to the stomach; the "food pipe."

### **Fistula; TE-fistula**

An abnormal duct or passage that connects an abscess or cavity to the body surface; a small opening between the trachea and the esophagus.

### **Indications**

A reason to prescribe a drug or perform a procedure.

### **Respiration**

The act or process of inhaling and exhaling; breathing; moving air into and out of the lungs. Also called ventilation.

---

**Respiratory Tract (laryngectomy)**

The air passage from the tracheostoma to the air sacs of the lungs, including the trachea, and bronchi through which air passes during breathing.

**Tracheoesophageal puncture (TEP)**

An opening made by a surgeon between the trachea and the esophagus. The opening is made at the time of initial surgery or later so that a voice prosthesis can be placed into this passage.

**Trachea**

The airway through which respiratory air passes to and from the lungs (windpipe).

**Tracheoesophageal (TE)**

Pertaining to the trachea and esophagus.

**Tracheostoma**

An opening through the neck and into the trachea. It is through this opening that the laryngectomye breathes.

**Voice prosthesis**

A one-way speech valve that is inserted into the TEP to allow air to pass from the trachea through the prosthesis and into the esophagus to produce sound.

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---

## BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFÍA / BIBLIOGRAFIA / BIBLIOGRAFIE

Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. / Se pueden solicitar más referencias. / Autres références disponibles sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. / Meer informatie is op aanvraag beschikbaar. / Referências adicionais mediante solicitação.

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.

- 
10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.
  11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
  12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
  13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
  14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
  15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
  16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

**INHEALTH**<sup>®</sup>  
T E C H N O L O G I E S

CE 0344

**InHealth Technologies**

A Division of Helix Medical, LLC  
1110 Mark Avenue  
Carpinteria, CA 93013-2918 USA  
[www.inhealth.com](http://www.inhealth.com)

**EMERGO EUROPE**

Molenstraat 15  
2513 BH, The Hague  
The Netherlands