INHEALTH™

TECHNOLOGIES

Blom-Singer® Indwelling TEP Occluder

USER Instructions For Use

Blom-Singer®

voice restoration systems

R1

37725-01E

37725-01E Effective January 2016 / Gültig ab Januar 2016 /
Efectivo a partir de enero de 2016 / Date d'entrée en vigueur :
janvier 2016 / Valido da gennaio 2016 / Met ingang van januari
2016 / Vigente em janeiro de 2016 / V platnosti od ledna 2016 /
Wchodzi w życie w styczniu 2016 r.

Blom-Singer and **InHealth Technologies** are registered
trademarks of **Freudenberg Medical, LLC**

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BLOM-SINGER® INDWELLING TEP OCCLUDER

Note: Some common terms your clinician will use are located in a glossary at the end of this chapter.


Diagram 1

Your clinician will place this device in your tracheostoma (diagram 1C) to keep your TE puncture tract open (diagram 1B). This device will also prevent liquids and food from entering your windpipe (diagram 1D) while eating, drinking, or swallowing. You will not be able to talk with this device.

This Blom-Singer medical device may be used only once by one patient. It may not be reused. See Warnings, Precautions and Complications.

Your clinician will talk to you about the risks and benefits of using this device. Your clinician will teach you how the device works. Your clinician will also show you how to care for your tracheostoma.

You should not use this instruction manual instead of talking with your clinician about the Indwelling TEP Occluder.

WHO SHOULD USE THIS DEVICE AND WHEN (INDICATIONS)

This device can only be used if your tracheostoma puncture is healed and undamaged. Only your clinician can place this device into your tracheostoma.

WHEN YOU SHOULD NOT USE THIS DEVICE (CONTRAINDICATIONS)

Your clinician will not be able to place this device into your tracheostoma if your stoma is too small or if the device makes it difficult for you to breathe. You will not be able to talk with this device.

HOW SUPPLIED

You will need a prescription from your clinician for this device (USA patients).

PRODUCT DESCRIPTION

This device is made of silicone, a material used in many medical products that is not harmful to your body. This device is latex-free.

WARNINGS AND PRECAUTIONS

Your clinician must provide to you instructions on how to use and care for this device. If there is something about this device that you do not completely understand, be sure to ask your clinician to explain.

Over time your health condition may change. You should have your clinician examine your device and TEP as soon as you notice any changes, for example:

- it has come out of your airway;
- it has tears or cracks;
- you are experiencing fluids leaking around the device;
- it doesn't fit the same in your TEP;
- it is not working properly.

You must have a TE device in your puncture **at all times**. Your TEP will close (sometimes within minutes or hours) if there is not a device in place to keep it open. If your TEP

closes, you may have to have an additional surgery.

If the Device Falls Out

If the Indwelling Occluder becomes loose from its original position (falls out), contact your clinician immediately. A foreign object in your airway can cause serious health issues such as an infection or pneumonia. Immediate signs of a loose device are gagging, choking, coughing, or trouble breathing.

If it is possible for you to locate where the device has become stuck in your body, do so immediately. If you think that the device has:

- become stuck in your windpipe, try to cough the Indwelling Occluder out. Immediately seek medical attention if you cannot get the device out of your windpipe by coughing;
- has fallen into your lungs, seek emergency medical attention immediately and contact your clinician;

- gone into your food pipe, you will likely pass it with a bowel movement. (If you have a stomach or intestinal health condition contact your clinician immediately if you note any changes, i.e., stomach pain.)

Using a mirror, insert a Blom-Singer puncture dilator, sold separately, or a suitable device as recommended by your clinician, into your puncture site. This will help to keep the puncture from closing and leaking fluids. Schedule an appointment with your clinician to insert a new device as soon as possible. A puncture dilator is available from your clinician with instructions on how to use it.

How Should You Prepare For A Planned or Unplanned Medical Situation?

Medical personnel will need to know about your special breathing needs so that they can give you proper care. Discuss with your clinician how you should identify yourself as a neck breather in the event of

an emergency or surgical procedure, and especially if you will need anesthesia. A medical alert bracelet or medical I.D. card is helpful and should be kept in a clearly visible location.

Device Insertion/Removal

Only a clinician, or qualified, trained medical professional, can insert and replace the Indwelling TEP Occluder for you.

The average lifetime of the Indwelling TEP Occluder varies from person to person. It also depends on how the device is used, and the conditions it is used under. Regular visits to your clinician are recommended every 6 months.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or a licensed speech pathologist.

COMPLICATIONS

Although rare, the following complications have been identified to occur with silicone

devices of the Blom-Singer type. Be sure to talk to your clinician about medical treatment if you experience any of these symptoms:

- pain, redness, heat, swelling, or infection around your tracheostoma or airway. These symptoms may require removal of the device and/or medical treatment;
- you notice a change in the status of your tissue, color/odor of sputum, and/or you have fever/flu-like symptoms, which could indicate the presence of an infection;
- accidental inhalation or swallowing of the device into your airway, food pipe (diagram 1A) and/or digestive tract, which may require removal by a doctor;
- fluids leaking around the device that may cause you to cough. If this condition does not go away you may have to have a surgery to revise your puncture. Your clinician may have to close your puncture;

-
- the device has become loose out of position in your TEP. This will require replacement of the device, and/or closure of your TE puncture;
 - the size of your TE puncture changes and becomes too big. This may cause fluids to leak around the device;
 - your body's response to a foreign object may cause irritated tissue to build-up around the puncture. This irritated tissue can make a "donut ring" or scarring around the prosthesis. These conditions may cause the device to rotate easily in the TEP;
 - difficulty swallowing;
 - itching, swelling, chest tightness, rash, nausea, vomiting, diarrhea or difficulty breathing following insertion of the device, seek medical treatment;
 - fever, low blood pressure, dizziness, rapid breathing, chills and shaking, or other flu like symptoms, seek medical treatment.

CLEANING AND CARE INSTRUCTIONS

The Indwelling TEP Occluder and user's hands **should always be thoroughly** clean to prevent contaminants from entering into the tracheostoma, puncture, or food pipe.

Note: It is not necessary to clean this device. Occasionally, on the portion of the device that shows out of your puncture you may need to remove dried mucus using tweezers. If you experience a lot of dried mucus, under the care of your clinician, you may want to consider use of humidification to soften the dried substance before removing. This may help to make the removal easier for you. Always be careful to avoid any liquids from entering your tracheostoma, as this will cause you to cough. **Do not attempt to use any liquids to clean the device.**

CONSUMER AFFAIRS

USA

If you have any questions or dissatisfaction with a product please contact our customer

service department via telephone, fax, post, or E-mail: TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594; POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, E-MAIL: productcomplaints@inhealth.com

International

Please contact your clinician or the InHealth distributor for your region.

GLOSSARY

Airway (Diagram 1D)

The air passage by which air reaches the lungs.

Clinician

A health professional, such as a physician, licensed speech-language pathologist, or nurse who is involved in clinical practice.

Contraindications

A reason that makes it inadvisable to

prescribe a particular procedure or treatment.

Esophagus (Diagram 1A)

The muscular membranous tube for the passage of food from the pharynx to the stomach; the “food pipe.”

Indications

A reason to prescribe a drug or perform a procedure.

Tracheoesophageal puncture (TEP) (Diagram 1B)

An opening made by a surgeon between the trachea and the esophagus. The opening is made at the time of initial surgery or later so that a prosthetic device can be placed into this passage.

Tracheostoma (Diagram 1C)

An opening through the neck and into the trachea. It is through this opening that the laryngectomee breathes.

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Additional references available upon request.

Weitere Referenzen aus der Literatur sind auf Anfrage erhältlich. / Se pueden solicitar más referencias. / Autres références disponibles sur demande. / Riferimenti aggiuntivi sono disponibili a richiesta. / Meer informatie is op aanvraag beschikbaar. / Referências adicionais mediante solicitação. / Další odkazy na literaturu jsou k dispozici na požádání. / Dodatkowe informacje dostępne na żądanie.

1. Blom, E.D.: "Some Comments on the Escalation of Tracheoesophageal Voice Prosthesis Dimensions," **Arch Otolaryngol**, 129: 500-502, 2003.
2. Eerenstein, S.E., Grolman, W., Schouwenburg, P.F.: "Downsizing of Voice Prosthesis Diameter in Patients with Laryngectomy," **Arch Otolaryngol**, 128: 838-841, 2003.
3. Hamaker, R.C., Blom, E.D.: "Botulinum Neurotoxin for Pharyngeal Constrictor Muscle Spasm in Tracheoesophageal Voice Restoration," **Laryngoscope**, 113:1479-1482, 2003.
4. Evans, P.H.R., Blom, E.D.: "Functional Restoration of Speech," **Principles and Practice of Head and Neck Oncology**, Evans, P.H.R., Montgomery, P.W., Gullane, P.J. (Eds.) Taylor and Francis Group, London, 2003.
5. Blom, E.D., Singer M.I., and Hamaker, R.C., **Tracheoesophageal Voice Restoration Following Total Laryngectomy**. San Diego: Singular Publishing Group, Inc.1998.
6. Blom, E.D., Hamaker, R.C.: "Tracheoesophageal Voice Restoration Following Total Laryngectomy," **Cancer of the Head and Neck**, E.N. and Suen, J. (Eds). W.B. Saunders Publishers. 839-852, 1996.
7. Blom, E.D., "Tracheoesophageal Speech," **Seminars in Speech and Language**. New York: Thieme Medical Publishers. 16:3, 191-204, 1995.
8. Pauloski, B.R., Blom, E.D., Logemann, J.A., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part II: Swallow Characteristics," **Laryngoscope**, 105:10: 1104-1110, 1995.
9. Blom, E.D., Pauloski, B.R., Hamaker, R.C. "Functional Outcome After Surgery for Prevention of Pharyngospasms in Tracheoesophageal Speakers. Part I: Speech Characteristics," **Laryngoscope**, 105:10: 1093-1103, 1995.
10. Blom, E.D., Hamaker, R.C., Freeman, S.B.: "Postlaryngectomy Voice Restoration," Lucente, F.E. (Ed.): **Highlights of the Instructional Courses**, Vol. 7. St. Louis: Mosby-Year Book. pp 3-10, 1994.

-
11. Blom, E.D., Singer M.I., and Hamaker, R.C. "A Prospective Study of Tracheoesophageal Speech," **Arch Otolaryngol**, 112: 440-447, 1986.
 12. Blom, E.D., Singer, M.I. "Disinfection of Silicone Voice Prostheses," **Arch Otolaryngol**, 112: 1986.
 13. Blom, E.D., Singer M.I., and Hamaker, R.C. "An Improved Esophageal Insufflation Test," **Arch Otolaryngol**, 111: 211-212, 1985.
 14. Hamaker, R.C., Singer M.I., Blom E.D. and Daniels, H.A. "Primary Voice Restoration at Laryngectomy," **Arch Otolaryngol**, 111: 182-186, 1985.
 15. Singer, M.I. and Blom, E.D. "A Selective Myotomy for Voice Restoration After Total Laryngectomy," **Arch Otolaryngol**, 107: 670-673, 1981.
 16. Singer, M.I. and Blom, E.D. "An Endoscopic Technique for Restoration of Voice After Laryngectomy," **Ann Otol Rhinol Laryngol**, 89:6: 529-533, 1980.

INHEALTH™

TECHNOLOGIES

UNITED STATES

1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
Toll-Free: 800.477.5969
Toll-Free Fax: 888.371.1530

INTERNATIONAL

Tel: +1.805.684.9337
Fax: +1.805.684.8594
E-mail: order@inhealth.com

ORDER ONLINE

www.inhealth.com

Freudenberg Medical, LLC
1110 Mark Avenue
Carpinteria, CA 93013-2918 USA
www.inhealth.com

EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands