

NOTICE

BLOM-SINGER HME SYSTEM USERS

TO AVOID DIMINISHED VOICE QUALITY AND AIR LEAKAGE USE ONLY WITH THE BLOM-SINGER ATTACHMENT SYSTEMS LISTED BELOW

| <i>Description</i> | <i>Catalog Number</i> |
|--|-----------------------|
| Blom-Singer TruSeal® Adhesive Housing, 7 pack | BE6053 |
| Blom-Singer TruSeal® Adhesive Housing, 30 pack | BE6054 |
| Tracheostoma Valve Housing, standard | BE6038 |
| Tracheostoma Valve Housing, large | BE6039 |
| Tracheostoma Valve Housing, standard, PVC | BE6038-P |
| Tracheostoma Valve Housing, large, PVC | BE6039-P |

THE BLOM-SINGER HME SYSTEM IS ALSO COMPATIBLE WITH THE BARTON-MAYO BUTTON

COMPATIBILITY WITH OTHER MANUFACTURERS ATTACHMENT SYSTEMS CANNOT BE ASSURED. DIMINISHED VOICE QUALITY AND AIR LEAKAGE MAY OCCUR.